



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SYNPO2 Polyclonal Antibody

多克隆抗体

产品货号: PA11210

第三版

描述: Synaptopodin-2 is a protein that in humans is encoded by the SYNPO2 gene. Has an actin-binding and actin-bundling activity. Synaptopodin is a proline-rich, actin-associated protein expressed in mature dendritic spines and renal podocytes. Synaptopodin appears to play a role in the actin-based plasticity of spines by linking actin to the spine apparatus.

抗原: Synthetic peptide of human SYNPO2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Skeletal muscle-specific.