

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse FOXD3 Polyclonal Antibody

多克隆抗体

产品货号: PA11275

第三版

描述:Forkhead box D3 also known as FOXD3 is a forkhead protein that in humans is encoded by the FOXD3 gene. This gene belongs to the forkhead protein family of transcription factors which is characterized by a DNA-binding forkhead domain. FoxD3 functions as a transcriptional repressor and contains the C-terminal engrailed homology-1 motif (eh1), which provides an interactive surface with a transcriptional co-repressor Grg4 (Groucho-related gene-4).

抗原:Synthetic peptide of human FOXD3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in chronic myeloid leukemia, Jurkat T-cell leukemia and teratocarcinoma cell lines, but not in any other cell lines or normal tissues examined.