

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GJA9 Polyclonal Antibody

多克隆抗体

产品货号: PA11276

第三版

描述:Connexins, such as GJA9, are involved in the formation of gap junctions, intercellular conduits that directly connect the cytoplasm of contacting cells. Each gap junction channel is formed by docking of 2 hemichannels, each of which contains 6 connexin subunits.

抗原:Synthetic peptide of human GJA9

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly abundant in skeletal muscle. Also detected in testis.