

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TGIF2LX Polyclonal Antibody

多克隆抗体

产品货号: PA11352

第三版

描述:This gene encodes a member of the TALE/TGIF homeobox family of transcription factors. Testis-specific expression suggests that this gene may play a role in spermatogenesis. A homolog of this gene lies within the male specific region of chromosome Y, in a block of sequence that is thought to be the result of a large X-to-Y transposition.

抗原:Recombinant protein of human TGIF2LX

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Specifically expressed in adult testis.