

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Bovine serum albumin BSA Polyclonal Antibody

多克隆抗体

产品货号: PA11613

第三版

描述: Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca²⁺, Na⁺, K⁺, fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.

抗原: Recombinant protein of human BSA

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Plasma.