

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ARMC10 Polyclonal Antibody

多克隆抗体

产品货号: PA11645

第三版

**描述:**This gene encodes a protein that contains an armadillo repeat and transmembrane domain. The encoded protein decreases the transcriptional activity of the tumor suppressor protein p53 through direct interaction with the DNA-binding domain of p53, and may play a role in cell growth and survival. Upregulation of this gene may play a role in hepatocellular carcinoma. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 3.

**抗原:**Recombinant protein of human ARMC10

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in all tissues tested with higher expression in placenta, liver, kidney, heart and brain.