

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SLC25A26 Polyclonal Antibody

多克隆抗体

产品货号: PA11714

第三版

**描述:** Mitochondrial carriers, including SLC25A26, are a family of transport proteins found mostly in the inner membranes of mitochondria. They shuttle metabolites and cofactors through the mitochondrial membrane

**抗原:** Recombinant protein of human SLC25A26

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Widely expressed. Highly expressed in testis, with moderate expression in brain, heart, kidney, lung, skeletal muscle, pancreas, small intestine and liver, and low expression in spleen.