



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CEACAM21 Polyclonal Antibody

多克隆抗体

产品货号: PA11728

第三版

**描述:**CEACAM21 is a gene which encodes carcinoembryonic antigen. It belongs to the immunoglobulin superfamily-CEA family. Three isoforms of the human protein are produced by alternative splicing. High levels of expression in hematopoietic system, brain and muscle.

**抗原:**Recombinant protein of human CEACAM21

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**