



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human S100PBP Polyclonal Antibody

多克隆抗体

产品货号: PA11742

第三版

描述: This gene encodes a protein that was originally identified by its interaction with S100 calcium-binding protein P. Expression of this protein has been reported to be associated with pancreatic ductal adenocarcinoma. Alternatively spliced transcript variants have been described.

抗原: Recombinant protein of human S100PBP

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in brain, spleen, and lung. Not detected in pancreas or liver. In pancreas, expressed predominantly in islet cells and to a lesser extent in acinar cells, but not expressed in ductal cells. Up-regulated in various pancreatic ductal adenocarcinomas and pancreatic intraepithelial neoplasias. Detected in pancreatic ductal adenocarcinoma cells (at protein level). Not detected in non-neoplastic ductal epithelium (at protein level).