

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRIM31 Polyclonal Antibody

多克隆抗体

产品货号: PA11763

第三版

描述:The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to both the cytoplasm and the nucleus. Its function has not been identified.

抗原:Recombinant protein of human TRIM31

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Up-regulated in gastric adenocarcinomas.1 Publication