
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RNASET2 Polyclonal Antibody

多克隆抗体

产品货号: PA11769

第三版

描述: This ribonuclease gene is a novel member of the Rh/T2/S-glycoprotein class of extracellular ribonucleases. It is a single copy gene that maps to 6q27, a region associated with human malignancies and chromosomal rearrangement.

抗原: Recombinant protein of human RNASET2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitous. Higher expression levels observed in the temporal lobe and fetal brain.