

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TMPRSS2 Polyclonal Antibody

多克隆抗体

产品货号: PA11805

第三版

描述: This gene encodes a protein that belongs to the serine protease family. The encoded protein contains a type II transmembrane domain, a receptor class A domain, a scavenger receptor cysteine-rich domain and a protease domain. Serine proteases are known to be involved in many physiological and pathological processes. This gene was demonstrated to be up-regulated by androgenic hormones in prostate cancer cells and down-regulated in androgen-independent prostate cancer tissue.

抗原: Recombinant protein of human TMPRSS2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed strongly in small intestine. Also expressed in prostate, colon, stomach and salivary gland.