

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GNL3 Polyclonal Antibody

多克隆抗体

产品货号: PA11848

第三版

描述:The protein encoded by this gene may interact with p53 and may be involved in tumorigenesis. The encoded protein also appears to be important for stem cell proliferation. This protein is found in both the nucleus and nucleolus. Three transcript variants encoding two different isoforms have been found for this gene.

抗原:Recombinant protein of human GNL3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Increased levels in lung tissue in cancer patients.