

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NRIP3 Polyclonal Antibody

多克隆抗体

产品货号: PA11851

第三版

描述:NRIP3 (nuclear receptor interacting protein 3), also known as C11orf14 or NY-SAR-105, is a 241 amino acid protein that is encoded by a gene which Human, Mouse maps to human chromosome 11. With approximately 135 million base pairs and 1,400 genes, chromosome 11 comprises approximately 4% of human genomic DNA and is considered a gene and disease association dense chromosome. The chromosome 11 encoded Atm gene is important for regulation of cell cycle arrest and apoptosis following double strand DNA breaks.

抗原:Recombinant protein of human NRIP3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: