

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NDUFA3 Polyclonal Antibody

多克隆抗体

产品货号: PA11871

第三版

描述:NDUFA3 (NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 3), also known as B9, is one of about 45 subunits comprising Complex I of the oxidative phosphorylation electron transport chain. Consisting of 84 amino acids and belonging to the Complex I NDUFA3 subunit family, NDUFA3 functions as an accessory subunit of the multi-protein mitochondrial membrane respiratory chain NADH dehydrogenase complex (known as Complex I). Complex I plays an important role in the transfer of electrons from NADH to the respiratory chain, a process that is essential for cellular respiration.

抗原:Recombinant protein of human NDUFA3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: