

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NDUFA9 Polyclonal Antibody

多克隆抗体

产品货号: PA11872

第三版

描述:The encoded protein is a subunit of the hydrophobic protein fraction of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. A pseudogene has been identified on chromosome 12.

抗原:Recombinant protein of human NDUFA9

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: