

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat BEX3 Polyclonal Antibody

多克隆抗体

产品货号: PA11882

第三版

**描述:**Protein BEX3 is a protein that in humans is encoded by the NGFRAP1 gene. May be a signaling adapter molecule involved in p75NTR-mediated apoptosis induced by NGF. Plays a role in zinc-triggered neuronal death By similarity. May play an important role in the pathogenesis of neurogenetic diseases.

**抗原:**Recombinant protein of human NGFRAP1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Found in ovarian granulosa cells, testis, prostate and seminal vesicle tissue. High levels also detected in liver.