

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MTFR1 Polyclonal Antibody

多克隆抗体

产品货号: PA11889

第三版

描述: This gene encodes a mitochondrial protein that is characterized by a poly-proline rich region. A chicken homolog of this protein promotes mitochondrial fission and the mouse homolog protects cells from oxidative stress. A related pseudogene of this gene is found on chromosome X.

抗原: Recombinant protein of human MTFR1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: