
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse MTF2 Polyclonal Antibody

多克隆抗体

产品货号: PA11890

第三版

描述:MTF2(MIER2(mesoderm induction early response 1, family member 2), also known as PCL2 or M96, is a protein that in humans is encoded by the MTF2 gene. It binds to the metal-regulating-element (MRE) of metallothionein gene promoter. May be involved in the activation of MT (metallothionein) genes in response to heavy-metal ions.

抗原:Recombinant protein of human MTF2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: