

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat MEP1A Polyclonal Antibody

多克隆抗体

产品货号: PA11911

第三版

**描述:**Meprin A subunit alpha also known as endopeptidase-2 or PABA peptide hydrolase is the alpha subunit of the meprin A enzyme that in humans is encoded by the MEP1A gene.[1][2]The MEP1A locus is on chromosome 6p in humans and on chromosome 17 in mice. The meprin alpha subunit product of the MEP1A gene is processed in the endoplasmic reticulum during intracellular transport, and is secreted as homomeric meprin A.

**抗原:**Recombinant protein of human MEP1A

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**