

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZC3H12A Polyclonal Antibody

多克隆抗体

产品货号: PA11917

第三版

描述:ZC3H12A is an MCP1 (CCL2; MIM 158105)-induced protein that acts as a transcriptional activator and causes cell death of cardiomyocytes, possibly via induction of genes associated with apoptosis.

抗原:Recombinant protein of human ZC3H12A

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in heart, placenta, spleen, kidney, liver and lung (PubMed:19909337). Expressed in leukocytes (PubMed:19909337). Expressed in monocyte (PubMed:16574901).