

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse KLHL1 Polyclonal Antibody

多克隆抗体

产品货号: PA11950

第三版

描述:KLHL1 (kelch-like protein 1) is a 748 amino acid protein that is the homolog of the Drosophila kelch protein. Localized to the cytoskeleton, KLHL1 contains six kelch repeats and one BTB (POZ) domain. The BTB (Broad-Complex, Tramtrack and Bric a brac) domain, also known as the POZ (Poxvirus and Zinc finger) domain, is an N-terminal homodimerization domain that contains multiple copies of kelch repeats and/or C2H2-type zinc fingers.

抗原:Recombinant protein of human KLHL1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: