

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human KCNK17 Polyclonal Antibody

多克隆抗体

产品货号: PA11960

第三版

**描述:**The protein encoded by this gene belongs to the family of potassium channel proteins containing two pore-forming P domains. This channel is an open rectifier which primarily passes outward current under physiological K<sup>+</sup> concentrations. This gene is activated at alkaline pH. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**抗原:**Recombinant protein of human KCNK17

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**