

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IRGM Polyclonal Antibody

多克隆抗体

产品货号: PA11977

第三版

描述: This gene encodes a member of the p47 immunity-related GTPase family. The encoded protein may play a role in the innate immune response by regulating autophagy formation in response to intracellular pathogens. Polymorphisms that affect the normal expression of this gene are associated with a susceptibility to Crohn's disease and tuberculosis.

抗原: Recombinant protein of human IRGM

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed (at protein level). Expressed in several tissues including colon, small bowel and peripheral blood leukocytes.