

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse TNFRSF12A Polyclonal Antibody

多克隆抗体

产品货号: PA12005

---

第三版

**描述:** Tumor necrosis factor receptor superfamily member 12A also known as the Tweak-receptor (TweakR) is a protein that in humans is encoded by the TNFRSF12A gene. Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins. Highly expressed in heart, placenta and kidney. Intermediate expression in lung, skeletal muscle and pancreas.

**抗原:** Recombinant protein of human TNFRSF12A

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Highly expressed in heart, placenta and kidney. Intermediate expression in lung, skeletal muscle and pancreas.