

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HSD17B2 Polyclonal Antibody

多克隆抗体

产品货号: PA12023

第三版

描述: Estradiol 17-beta-dehydrogenase 2 is an enzyme that in humans is encoded by the HSD17B2 gene. Capable of catalyzing the interconversion of testosterone and androstenedione, as well as estradiol and estrone. Also has 20-alpha-HSD activity. Uses NADH while EDH17B3 uses NADPH.

抗原: Recombinant protein of human HSD17B2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: