

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat HSD17B7 Polyclonal Antibody

多克隆抗体

产品货号: PA12024

---

第三版

**描述:** HSD17B7 encodes an enzyme that functions both as a 17-beta-hydroxysteroid dehydrogenase in the biosynthesis of sex steroids and as a 3-ketosteroid reductase in the biosynthesis of cholesterol.

**抗原:** Recombinant protein of human HSD17B7

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Highly expressed in adrenal gland, liver, lung and thymus. Expressed in breast, ovaries, pituitary gland, pregnant uterus, prostate, kidney, lymph node, small intestine, spinal cord and trachea. Weakly expressed in all other tissues tested. Isoform 3 is expressed in eye ciliary epithelial cells and neuroendocrine cells.