

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse GABPB2 Polyclonal Antibody

多克隆抗体

产品货号: PA12075

第三版

描述:GA-binding protein subunit beta-2 is a protein that in humans is encoded by the GABPB2 gene. The GABPB variant containing a C-terminal leucine zipper motif, which they called GABP-beta, can homodimerize and form a transcriptionally active heterotetrameric complex with GABPA. The variant lacking this motif, which they called GABP-gamma, can form heterodimers but cannot form transcriptionally active heterotetramers with GABPA.

抗原:Recombinant protein of human GABPB2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: