

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FAAH2 Polyclonal Antibody

多克隆抗体

产品货号: PA12122

第三版

描述: This gene encodes a fatty acid amide hydrolase that shares a conserved protein motif with the amidase signature family of enzymes. The encoded enzyme is able to catalyze the hydrolysis of a broad range of bioactive lipids, including those from the three main classes of fatty acid amides; N-acyl ethanolamines, fatty acid primary amides and N-acyl amino acids. This enzyme has a preference for monounsaturated acyl chains as a substrate.

抗原: Recombinant protein of human FAAH2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in the brain, small intestine and testis. Also expressed in the heart, kidney, liver, lung and prostate.