

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EPST11 Polyclonal Antibody

多克隆抗体

产品货号: PA12139

第三版

**描述:**EPST11 (Epithelial-stromal interaction protein 1), also known as BRES11, is a 318 amino acid protein that is highly expressed in spleen, placenta, small intestine, liver, kidney, thymus, testis and salivary gland. Existing as three isoforms that are formed due to alternative splicing events, EPST11 is weakly expressed in normal breast tissue, but is heavily up-regulated in breast carcinoma, suggesting an important role for EPST11 in tumor formation and/or progression. Additionally, EPST11 is found in blood cells from systemic lupus erythematosus (SLE)-afflicted patients, implicating EPST11 as a potential protein involved in SLE.

**抗原:**Recombinant protein of human EPST11

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**