

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EMD Polyclonal Antibody

多克隆抗体

产品货号: PA12145

第三版

**描述:**Emerin is a serine-rich nuclear membrane protein and a member of the nuclear lamina-associated protein family. It mediates membrane anchorage to the cytoskeleton. Drefuss-Emery muscular dystrophy is an X-linked inherited degenerative myopathy resulting from mutation in the emerin gene.

**抗原:**Recombinant protein of human EMD

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Skeletal muscle, heart, colon, testis, ovary and pancreas.