

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ECE2 Polyclonal Antibody

多克隆抗体

产品货号: PA12146

第三版

**描述:**Endothelin-converting enzymes, such as ECE2 (EC 3.4.24.71), are type II metalloproteases that generate functionally pleiotropic members of the endothelin vasoactive peptide family. It converts big endothelin-1 to endothelin-1. Also involved in the processing of various neuroendocrine peptides, including neurotensin, angiotensin I, substance P, proenkephalin-derived peptides, and prodynorphin-derived peptides. May limit beta-amyloid peptide accumulation in brain. May also have methyltransferase activity.

**抗原:**Recombinant protein of human ECE2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**Expressed in brain. Strongly down-regulated in inferior parietal lobe from Alzheimer disease patients (at protein level).