

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FOXR2 Polyclonal Antibody

多克隆抗体

产品货号: PA12202

第三版

描述:Forkhead-box (FOX) transcription factors are implicated in carcinogenesis through gene amplification, retroviral integration, or chromosomal translocation. It expressed in breast cancer cell lines and primary cancer. The Subcellular location is in nucleus.

抗原:Recombinant protein of human FOXR2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in breast cancer cell lines and primary cancer.