

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat DDIT4L Polyclonal Antibody

多克隆抗体

产品货号: PA12238

第三版

描述:REDD-2 (regulated in development and DNA damage response 2), also designated Rtp801L or DDIT4L (DNA-damage-inducible transcript 4-like), is a 193 amino acid cytoplasmic protein belonging to the DDIT4 family and is predominantly expressed in skeletal muscle. Considered a stress-induced protein, REDD-2 is a negative regulator of the mTOR (mammalian target of rapamycin) pathway. mTOR is a serine/threonine kinase that plays an essential role in cell growth control and is an important regulator of skeletal muscle size. Highly expressed in human atherosclerotic lesions and macrophages, REDD-2 mediates monocyte cell death through reduction of Trx (thioredoxin-1) expression. REDD2 expression in macrophages increases oxidized LDL (oxLDL)-induced cell death, suggesting that REDD2 may play a critical role in arterial pathology.

抗原:Recombinant protein of human DDIT4L

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Up-regulated in atherosclerotic plaques relative to healthy segments of the same artery.