

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CTAGE5 Polyclonal Antibody

多克隆抗体

产品货号: PA12310

第三版

**描述:**The protein encoded by this gene is a tumor-associated antigen found in cutaneous T-cell lymphoma and several other cancers. Autoantibodies against the encoded protein have been found in some cancers. Several transcript variants encoding different isoforms have been found for this gene.

**抗原:**Recombinant protein of human CTAGE5

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Isoform 5A is expressed only in testis (at the protein level). Other isoforms are expressed in several other normal tissues, including brain, muscle and cranial skin. Isoform 5A (at protein level) and isoform 5B are expressed in cutaneous T-cell lymphoma (CTCL) cell lines, colorectal carcinomas, breast carcinomas and melanoma. Isoform 5B, but not isoform 5A, is expressed in head and neck squamous cell carcinoma.