

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IFNG Polyclonal Antibody

多克隆抗体

产品货号: PA12325

第三版

描述: This gene encodes a member of the type II interferon family. The protein encoded is a soluble cytokine with antiviral, immunoregulatory and anti-tumor properties and is a potent activator of macrophages. Mutations in this gene are associated with aplastic anemia.

抗原: Recombinant protein of human IFNG

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Released primarily from activated T lymphocytes.