

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CKAP2 Polyclonal Antibody

多克隆抗体

产品货号: PA12351

第三版

描述: This gene encodes a cytoskeleton-associated protein that stabilizes microtubules and plays a role in the regulation of cell division. The encoded protein is itself regulated through phosphorylation at multiple serine and threonine residues. There is a pseudogene of this gene on chromosome 14. Alternative splicing results in multiple transcript variations.

抗原: Recombinant protein of human CKAP2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Abundant in testis, thymus, and in tumor derived cell lines, while barely detectable in liver, prostate, and kidney.