

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NOS2 Polyclonal Antibody

多克隆抗体

产品货号: PA12509

第三版

**描述:** Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. This gene encodes a nitric oxide synthase which is expressed in liver and is inducible by a combination of lipopolysaccharide and certain cytokines. Three related pseudogenes are located within the Smith-Magenis syndrome region on chromosome 17.

**抗原:** Recombinant protein of human NOS2

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expressed in the liver, retina, bone cells and airway epithelial cells of the lung. Not expressed in the platelets.