

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLAMF8 Polyclonal Antibody

多克隆抗体

产品货号: PA12519

第三版

**描述:**This gene encodes a member of the CD2 family of cell surface proteins involved in lymphocyte activation. These proteins are characterized by Ig domains. This protein is expressed in lymphoid tissues, and studies of a similar protein in mouse suggest that it may function during B cell lineage commitment. The gene is found in a region of chromosome 1 containing many CD2 genes.

**抗原:**Recombinant protein of human SLAMF8

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in lymph node, spleen, thymus and bone marrow.