

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human KCNMA1 Polyclonal Antibody

多克隆抗体

产品货号: PA12523

第三版

描述:MaxiK channels are large conductance, voltage and calcium-sensitive potassium channels which are fundamental to the control of smooth muscle tone and neuronal excitability. MaxiK channels can be formed by 2 subunits: the pore-forming alpha subunit, which is the product of this gene, and the modulatory beta subunit. Intracellular calcium regulates the physical association between the alpha and beta subunits. Alternatively spliced transcript variants encoding different isoforms have been identified.

抗原:Recombinant protein of human KCNMA1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed. Except in myocytes, it is almost ubiquitously expressed.