

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PIWIL4 Polyclonal Antibody

多克隆抗体

产品货号: PA12590

第三版

描述:PIWIL4 belongs to the Argonaute family of proteins, which function in development and maintenance of germline stem cells. It plays a central role during spermatogenesis by repressing transposable elements and prevent their mobilization, which is essential for the germline integrity. May be involved in the chromatin-modifying pathway by inducing 'Lys-9' methylation of histone H3 at some loci.

抗原:Recombinant protein of human PIWIL4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in testis. According to PubMed:17544373, it is ubiquitously expressed.