

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PRL Polyclonal Antibody

多克隆抗体

产品货号: PA12626

第三版

描述:This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein.

抗原:Recombinant protein of human PRL

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: