

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NOL9 Polyclonal Antibody

多克隆抗体

产品货号: PA12655

第三版

**描述:** Nucleolar protein 9 is a protein that in humans is encoded by the NOL9 gene. Polynucleotide 5'-kinase involved in rRNA processing. The kinase activity is required for the processing of the 32S precursor into 5.8S and 28S rRNAs, more specifically for the generation of the major 5.8S(S) form. In vitro, has both DNA and RNA 5'-kinase activities. Probably binds RNA.

**抗原:** Recombinant protein of human NOL9

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**