

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SFTPD Polyclonal Antibody

多克隆抗体

产品货号: PA12716

第三版

**描述:** Surfactant, pulmonary-associated protein D, also known as SFTPD or SP-D, is a protein which in humans is encoded by the SFTPD gene. SFTPD is an innate immune system collectin. Surfactant protein D has been shown to interact with DMBT1. Contributes to the lung's defense against inhaled microorganisms. May participate in the extracellular reorganization or turnover of pulmonary surfactant. Binds strongly maltose residues and to a lesser extent other alpha-glucosyl moieties.

**抗原:** Recombinant protein of human SFTPD

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expressed in lung, brain, pancreas and adipose tissue (mainly mature adipocytes).