

---

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat HIPK4 Polyclonal Antibody

多克隆抗体

产品货号: PA12734

---

第三版

**描述:** HIPK4 (Homeodomain-interacting protein kinase 4) is a 616 amino acid protein that localizes to the cytoplasm and exists as two isoforms that are produced by alternative splicing events. Like other members of the HIPK family, HIPK4 functions as a protein kinase that catalyzes the ATP-dependent phosphorylation of target proteins and is thought to act as a corepressor of target transcription factors.

**抗原:** Recombinant protein of human HIPK4

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**