

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SYT7 Polyclonal Antibody

多克隆抗体

产品货号: PA12735

第三版

描述: This gene is a member of the synaptotagmin gene family and encodes a protein similar to other family members that mediate calcium-dependent regulation of membrane trafficking in synaptic transmission. A similar protein in rodents mediates hormone secretion and lysosome exocytosis. In humans, expression of this gene has been associated with prostate cancer. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

抗原: Recombinant protein of human SYT7

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in a variety of adult and fetal tissues.