

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SPATA7 Polyclonal Antibody

多克隆抗体

产品货号: PA12741

第三版

**描述:** This gene, originally isolated from testis, is also expressed in retina. Mutations in this gene are associated with Leber congenital amaurosis and juvenile retinitis pigmentosa. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**抗原:** Recombinant protein of human SPATA7

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**