

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SULT1E1 Polyclonal Antibody

多克隆抗体

产品货号: PA12745

第三版

**描述:** Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. The gene structure (number and length of exons) is similar among family members. This gene encodes a protein that transfers a sulfo moiety to and from estrone, which may control levels of estrogen receptors.

**抗原:** Recombinant protein of human SULT1E1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Liver, intestine and at lower level in the kidney.