

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RPL26L1 Polyclonal Antibody

多克隆抗体

产品货号: PA12787

第三版

**描述:** This gene encodes a protein that shares high sequence similarity with ribosomal protein L26. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. Transcript variants utilizing alternative polyA signals exist.

**抗原:** Recombinant protein of human RPL26L1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:**