

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BCL2L14 Polyclonal Antibody

多克隆抗体

产品货号: PA12807

第三版

描述:The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene.

抗原:Recombinant protein of human BCL2L14

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Isoform 1 is widely expressed. Isoform 2 is testis-specific.